



**NEPHROLOGY PROGRAM
DEPARTMENT POLICIES AND PROCEDURES**

**Biomed Neph - Section 01- Equipment maintenance - Neph Tech 1-13
PM Program for the DWA CeCon 3000 Concentrate Supply Unit**

No.: 01468 (TOH Standardized Policy Number)

ISSUED BY: Nephrology Technical Practice Committee
DATE OF APPROVAL: 2015/03

APPROVED BY: Program Clinical Director / Division Head
LAST REVIEW/REVISION DATE: 2017/03

CATEGORY: Equipment Maintenance
IMPLEMENTATION DATE: 2015/03

POLICY STATEMENT:

- In order to ensure patient safety all DWA Cecon concentrate supply units shall go through a preventive maintenance program bi-yearly (every 6 months) following manufacturer specifications

DEFINITION(S): N/A

ALERTS: N/A

TOOLS/SUPPLIES REQUIRED:

- Safety glasses
- Standard hand tools
- 100ml beaker
- Pocket timer
- Multi-meter
- IBP or Mesa dialysate meter with pressure testing tool (pre made jig)

PROCEDURE:

1. The first technologist to initiate the preventive maintenance shall open a NephroCare work order as well as initiate the appropriate preventive maintenance checklist
2. All interventions related to this specific preventive maintenance shall be noted on the preventive maintenance checklist

3. Follow the DWA CeCon preventive maintenance instructions: Refer to Appendix A below
4. Complete the DWA CeCon preventive maintenance check list: Refer to Appendix B below
5. Electrical safety testing procedures must be carried out following Policy # 01061-Neph Tech 1-08 Electrical Safety Testing Procedure
6. ESA testing will only include protective earth resistance tests, and test values are indicated on the CeCon preventive maintenance check list
7. The technologist completing the maintenance shall be responsible to enter the final information in the NephroCare work order and closing it
8. The technologist completing the maintenance shall update the Equipment Maintenance Log Sheet on the "V" drive
9. The preventive maintenance checklist shall be placed in the appropriate equipment binder along with a copy of the NephroCare work order

RELATED POLICIES / LEGISLATION:

1. Nephrology Policies and Procedures - [Biomed Neph - Section 01 - Equipment Maintenance - Neph Tech 1-07 \(#01060\) NephroCare Work Orders](#)
2. Nephrology Policies and Procedures - [Biomed Neph - Section 01 - Equipment Maintenance - Neph Tech 1-08 \(#01061\) Electrical Safety Testing Procedure](#)

REFERENCE:

1. DWA gmbh & Co.KG Concentrate Supply Unit CeCon 3000 service manual (CASK1E-0 Maintenance-SC CeCon 3000)



DWA CeCon 3000 Functional Testing Instructions for Bi-Yearly PM Checklist (March 21, 2017)

- Remove the 6 screws from the top panel to gain access to the ST-3000 board. To secure the top panel, align the top middle hole of the panel with the front center hole of the Cecon housing and fasten it with one short screw. See Diagram 1
- For all calibration adjustments and electronic component locations see Diagram 1

Pressure Calibration

There are two differences between QCH and AKU.


- QCH pumps run only on demand and the check valves are located on the output of the pumps.
- AKU pumps are continuous feed and the check valves are located in the white distribution block on the wall.

See Diagram 2 for QCH set up


See Diagram 3 for AKU set up

1. Push the service button once (diag1)
2. using the control panel keys (diag 4) record the pressure set point values. PM Sec 4.0
Normally QCH = 0.5 bar (375mmHg)
AKU = 0.28 bar (210mmHg)
3. Connect the calibration tool as shown in diagram 2 and 3. Be sure to have transducer protector on the meter being used.
 - QCH – disconnect the sensor supply line from the input to the white distribution block (Diag.2)
 - AKU – disconnect the sensor return/drain line from the white sensor block (Diag.3)
4. Turn on system if not already on
5. Create a flow restriction using the shut off valve on the calibration tool.
 - Using a small beaker and timer create a flow. (fast drip).
 - QCH 30-40ml / min
 - AKU 100ml / min
6. Record the pressure for each line sensor. Pressure should be +/- 75 mmHg from set point. If necessary calibrate using Pressure sensor pots (diag1). Pots are indicated on the PM Checklist. (The pots are very sensitive adjustments)

Diag.1




Diag.2 (QCH)




This is a controlled document prepared solely for use at The Ottawa Hospital (TOH). TOH accepts no responsibility for use of this material by any person or organization not associated with TOH. No part of this document may be reproduced in any form for publication without permission of TOH. A printed copy may not reflect the current electronic document and should be checked against the one on the TOH Intranet.

Diag.3



Diag.4



This is a controlled document prepared solely for use at The Ottawa Hospital (TOH). TOH accepts no responsibility for use of this material by any person or organization not associated with TOH. No part of this document may be reproduced in any form for publication without permission of TOH. A printed copy may not reflect the current electronic document and should be checked against the one on the TOH Intranet.

4. Operating Values and Pressure Verification:

4.1 Parameters

Factor: P _____ time vent: _____ sec
I _____ Time vent. A1: _____ sec
D _____ Time vent. A2: _____ min
N _____ TA: _____ sec (QCH)

Compare all parameters from previous PM. (these values should not change)..... _____ (initial)**

4.2 Pressure and Pump

QCH-0.5bar=375mmHg AKU 0.28bar=210mmhg

	Press1	Pres 2	Pres3	Pres4	
Set point:	_____	_____	_____	_____	Bar
Measured pres:	_____	_____	_____	_____	mmHg +/-75 from set point
Calibrated pres:	_____	_____	_____	_____	mmHg mark N/A if no adjustment
	P5	P6	P7	P8	
Pump Capacity after Pres.calibration:	_____	_____	_____	_____	%
Pump Capacity last PM:	_____	_____	_____	_____	%

4.3 Control Parameters: QCH only

Press.Diff: _____
Pump Tol: _____
Factor Na: _____

5. Clean CeCon: All internal, external, tubing and fittings..... _____ (initial)